

Product datasheet for GA101305

Human Dishevelled 3 (DVL3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DVL3 CRISPRa kit - CRISPR gene activation of human dishevelled segment polarity protein 3
Symbol:	DVL3
UniProt ID:	Q92997
Synonyms:	DRS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)