

Product datasheet for **GA101326**

Human EDG3 (S1PR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S1PR3 CRISPRa kit - CRISPR gene activation of human sphingosine-1-phosphate receptor 3
Symbol:	S1PR3
UniProt ID:	Q99500
Synonyms:	EDG-3; EDG3; LPB3; S1P3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)