

Product datasheet for **GA101350**

Human Flamingo homolog 1 (CELSR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CELSR3 CRISPRa kit - CRISPR gene activation of human cadherin EGF LAG seven-pass G-type receptor 3
Symbol:	CELSR3
UniProt ID:	Q9NYQ7
Synonyms:	ADGRC3; CDHF11; EGFL1; FMI1; HFMI1; MEGF2; RESDA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)