

## Product datasheet for **GA101361**

### Human eIF1A (EIF1AX) activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | EIF1AX CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 1A X-linked |
| Symbol:       | EIF1AX  |
| UniProt ID:   | <u><a href="#">P47813</a></u>   |
| Synonyms:     | eIF-1A; eIF-4C; EIF1A; EIF1AP1; EIF4C   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                   |

#### Product images:



[View online >](#)