


Product datasheet for **GA101364**

Human EIF2G (EIF2S3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF2S3 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 2 subunit gamma
Symbol:	EIF2S3
Components:	GA101364G1, EIF2G gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101364G2, EIF2G gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101364G3, EIF2G gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P41091</u>
Synonyms:	eIF-2gA; EIF2; EIF2G; EIF2gamma
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)