

Product datasheet for **GA101366**

Human EIF4A1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF4A1 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4A1
Symbol:	EIF4A1
UniProt ID:	P60842
Synonyms:	DDX2A; EIF-4A; eIF-4A-I; EIF4A; eIF4A-I
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)