

Product datasheet for **GA101367**

Human EIF4A2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF4A2 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4A2
Symbol:	EIF4A2
UniProt ID:	Q14240
Synonyms:	BM-010; DDX2B; eIF-4A-II; EIF4A; eIF4A-II; EIF4F
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)