

Product datasheet for **GA101395**

Human EMP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EMP1 CRISPRa kit - CRISPR gene activation of human epithelial membrane protein 1
Symbol:	EMP1
Components:	GA101395G1, EMP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101395G2, EMP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101395G3, EMP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P54849</u>
Synonyms:	CL-20; EMP-1; TMP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)