

## Product datasheet for **GA101412**

### Human EP300 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	EP300 CRISPRa kit - CRISPR gene activation of human E1A binding protein p300
Symbol:	EP300
UniProt ID:	<a href="https://www.uniprot.org/entry/Q09472">Q09472</a>
Synonyms:	KAT3B; p300; RSTS2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)