

Product datasheet for **GA101456**

Human ERV31 (ERV3-1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ERV3-1 CRISPRa kit - CRISPR gene activation of human endogenous retrovirus group 3 member 1, envelope
Symbol:	ERV3-1
Components:	GA101456G1, ERV31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101456G2, ERV31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101456G3, ERV31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q14264
Synonyms:	envR; ERV-R; ERV3; ERVR; HERV-R; HERVR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)