

Product datasheet for **GA101474**

Human ETV3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | ETV3 CRISPRa kit - CRISPR gene activation of human ETS variant 3 |
| Symbol: | ETV3 |
| UniProt ID: | P41162 |
| Synonyms: | bA110J1.4; METS; PE-1; PE1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)