

Product datasheet for **GA101496**

Human Thrombin Receptor (F2R) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	F2R CRISPRa kit - CRISPR gene activation of human coagulation factor II thrombin receptor
Symbol:	F2R
UniProt ID:	P25116
Synonyms:	CF2R; HTR; PAR-1; PAR1; TR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)