

Product datasheet for **GA101497**

Human PAR2 (F2RL1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | F2RL1 CRISPRa kit - CRISPR gene activation of human F2R like trypsin receptor 1 |
| Symbol: | F2RL1 |
| UniProt ID: | <u>P55085</u> |
| Synonyms: | GPR11; PAR2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)