

## Product datasheet for **GA101498**

### Human Proteinase Activated Receptor 3 (F2RL2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	F2RL2 CRISPRa kit - CRISPR gene activation of human coagulation factor II thrombin receptor like 2
Symbol:	F2RL2
UniProt ID:	<u><a href="#">O00254</a></u>
Synonyms:	PAR-3; PAR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)