

Product datasheet for **GA101567**

Human Fibrinogen beta chain (FGB) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FGB CRISPRa kit - CRISPR gene activation of human fibrinogen beta chain
Symbol:	FGB
UniProt ID:	P02675
Synonyms:	HEL-S-78p
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)