


Product datasheet for **GA101600**

Human FKBP25 (FKBP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FKBP3 CRISPRa kit - CRISPR gene activation of human FKBP prolyl isomerase 3
Symbol:	FKBP3
UniProt ID:	Q00688
Synonyms:	FKBP-3; FKBP-25; FKBP25; PPIase
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)