

Product datasheet for **GA101654**

Human FPR3 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | FPR3 CRISPRa kit - CRISPR gene activation of human formyl peptide receptor 3 |
| Symbol: | FPR3 |
| UniProt ID: | P25089 |
| Synonyms: | FML2_HUMAN; FMLP-R-II; FMLPY; FPRH1; FPRH2; FPRL2; RMLP-R-I |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)