

Product datasheet for **GA101802**

Human Gastric Inhibitory Polypeptide Receptor (GIPR) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GIPR CRISPRa kit - CRISPR gene activation of human gastric inhibitory polypeptide receptor
Symbol:	GIPR
UniProt ID:	<u>P48546</u>
Synonyms:	PGQTL2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)