

Product datasheet for **GA101857**

Human GNG3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GNG3 CRISPRa kit - CRISPR gene activation of human G protein subunit gamma 3
Symbol:	GNG3
UniProt ID:	P63215
Synonyms:	guanine nucleotide-binding protein gamma-3 subunit; guanine nucleotide binding protein (G protein), gamma 3; NBP gamma-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)