

Product datasheet for **GA101861**

Human GNG10 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | GNG10 CRISPRa kit - CRISPR gene activation of human G protein subunit gamma 10 |
| Symbol: | GNG10 |
| UniProt ID: | P50151 |
| Synonyms: | guanine nucleotide binding protein (G protein), gamma 10; guanine nucleotide binding protein 10; OTTHUMP00000021914; OTTHUMP00000063926 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)