

Product datasheet for **GA101893**

Human GPR4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR4 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 4
Symbol:	GPR4
UniProt ID:	P46093
Synonyms:	G protein-coupled receptor 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)