

Product datasheet for **GA101898**

Human CXCR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CXCR3 CRISPRa kit - CRISPR gene activation of human C-X-C motif chemokine receptor 3
Symbol:	CXCR3
UniProt ID:	P49682
Synonyms:	CD182; CD183; CKR-L2; CMKAR3; GPR9; IP10-R; Mig-R; MigR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)