

Product datasheet for **GA101909**

Human P2Y9 (LPAR4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LPAR4 CRISPRa kit - CRISPR gene activation of human lysophosphatidic acid receptor 4
Symbol:	LPAR4
UniProt ID:	Q99677
Synonyms:	GPR23; LPA4; P2RY9; P2Y5-LIKE; P2Y9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)