

Product datasheet for **GA101917**

Human GPR33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR33 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 33 (gene/pseudogene)
Symbol:	GPR33
Components:	GA101917G1, GPR33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA101917G2, GPR33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA101917G3, GPR33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q49SQ1</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)