

Product datasheet for **GA101921**

Human GPR37 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | GPR37 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 37 |
| Symbol: | GPR37 |
| UniProt ID: | <u>Q15354</u> |
| Synonyms: | EDNRBL; hET(B)R-LP; PAELR |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)