

Product datasheet for **GA102043**

Human H10 (H1F0) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	H1F0 CRISPRa kit - CRISPR gene activation of human H1 histone family member 0
Symbol:	H1-0
Components:	GA102043G1, H10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102043G2, H10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102043G3, H10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P07305</u>
Synonyms:	H1FV; H10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

