

Product datasheet for **GA102045**

Human Histone H1.3 (HIST1H1D) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H1D CRISPRa kit - CRISPR gene activation of human histone cluster 1 H1 family member d
Symbol:	H1-3
UniProt ID:	P16402
Synonyms:	H1.3; H1D; H1F3; H1s-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)