


Product datasheet for **GA102046**

Human Histone H1.4 (HIST1H1E) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H1E CRISPRa kit - CRISPR gene activation of human histone cluster 1 H1 family member e
Symbol:	H1-4
Components:	GA102046G1, Histone H1.4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102046G2, Histone H1.4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102046G3, Histone H1.4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P10412</u>
Synonyms:	dj221C16.5; H1.4; H1E; H1F4; H1s-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)