

Product datasheet for **GA102047**

Human Histone H1.5 (HIST1H1B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H1B CRISPRa kit - CRISPR gene activation of human histone cluster 1 H1 family member b
Symbol:	H1-5
UniProt ID:	<u>P16401</u>
Synonyms:	H1; H1.5; H1B; H1F5; H1s-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)