

Product datasheet for **GA102048**

Human Histone H1t (HIST1H1T) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H1T CRISPRa kit - CRISPR gene activation of human histone cluster 1 H1 family member t
Symbol:	H1-6
UniProt ID:	P22492
Synonyms:	dj221C16.2; H1.6; H1FT; H1t
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)