

## Product datasheet for **GA102055**

### Human H3.3A (H3F3A) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	H3F3A CRISPRa kit - CRISPR gene activation of human H3 histone family member 3A
Symbol:	H3-3A
UniProt ID:	<a href="#">P84243</a>
Synonyms:	H3.3A; H3F3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)