

## Product datasheet for **GA102057**

### Human Histone H1.1 (HIST1H1A) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H1A CRISPRa kit - CRISPR gene activation of human histone cluster 1 H1 family member a
Symbol:	H1-1
Components:	GA102057G1, Histone H1.1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102057G2, Histone H1.1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102057G3, Histone H1.1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q02539</a></u>
Synonyms:	H1.1; H1A; H1F1; HIST1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)