

## Product datasheet for **GA102068**

### Human HAS3 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | HAS3 CRISPRa kit - CRISPR gene activation of human hyaluronan synthase 3 |
| Symbol:       | HAS3   |
| UniProt ID:   | <u><a href="#">O00219</a></u>  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online >](#)