

Product datasheet for **GA102112**

Human HIVEP2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | HIVEP2 CRISPRa kit - CRISPR gene activation of human human immunodeficiency virus type I enhancer binding protein 2 |
| Symbol: | HIVEP2 |
| UniProt ID: | P31629 |
| Synonyms: | HIV-EP2; MBP-2; MIBP1; SHN2; ZAS2; ZNF40B |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)