

Product datasheet for **GA102147**

Human HMG4 (HMGB3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HMGB3 CRISPRa kit - CRISPR gene activation of human high mobility group box 3
Symbol:	HMGB3
UniProt ID:	O15347
Synonyms:	HMG-2a; HMG-4; HMG2A; HMG4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)