

Product datasheet for **GA102155**

Human CD168 (HMMR) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | HMMR CRISPRa kit - CRISPR gene activation of human hyaluronan mediated motility receptor |
| Symbol: | HMMR |
| UniProt ID: | O75330 |
| Synonyms: | CD168; IHABP; RHAMM |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)