

## Product datasheet for **GA102177**

### Human HNRPH3 (HNRNPH3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HNRNPH3 CRISPRa kit - CRISPR gene activation of human heterogeneous nuclear ribonucleoprotein H3
Symbol:	HNRNPH3
UniProt ID:	<a href="#">P31942</a>
Synonyms:	2H9; HNRPH3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)