

Product datasheet for **GA102180**

Human hnRNP U (HNRNPU) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HNRNPU CRISPRa kit - CRISPR gene activation of human heterogeneous nuclear ribonucleoprotein U
Symbol:	HNRNPU
UniProt ID:	Q00839
Synonyms:	hnRNP U; HNRNPU-AS1; HNRPU; SAF-A; SAFA; U21.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)