

Product datasheet for **GA102214**

Human HOXD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HOXD3 CRISPRa kit – CRISPR gene activation of human homeobox D3
Symbol:	HOXD3
Components:	GA102214G1, HOXD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102214G2, HOXD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102214G3, HOXD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P31249</u>
Synonyms:	Hox-4.1; HOXD; HOX4; HOX4A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

