

Product datasheet for **GA102268**

Human IGSF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IGSF3 CRISPRa kit - CRISPR gene activation of human immunoglobulin superfamily member 3
Symbol:	IGSF3
UniProt ID:	<u>O75054</u>
Synonyms:	EWI-3; LCDD; V8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)