

## Product datasheet for **GA102278**

### Human HTN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HTN3 CRISPRa kit - CRISPR gene activation of human histatin 3
Symbol:	HTN3
UniProt ID:	<a href="#">P15516</a>
Synonyms:	HIS2; HTN2; HTN5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)