

Product datasheet for **GA102287**

Human 5HT3A receptor (HTR3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HTR3A CRISPRa kit - CRISPR gene activation of human 5-hydroxytryptamine receptor 3A
Symbol:	HTR3A
UniProt ID:	<u>P46098</u>
Synonyms:	5-HT-3; 5-HT3A; 5-HT3R; 5HT3R; HTR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)