

Product datasheet for **GA102293**

Human Tenascin C (TNC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNC CRISPRa kit - CRISPR gene activation of human tenascin C
Symbol:	TNC
UniProt ID:	P24821
Synonyms:	150-225; DFNA56; GMEM; GP; HXB; JI; TN; TN-C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)