

Product datasheet for **GA102307**

Human ID3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ID3 CRISPRa kit - CRISPR gene activation of human inhibitor of DNA binding 3, HLH protein
Symbol:	ID3
UniProt ID:	Q02535
Synonyms:	bHLHb25; HEIR-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)