

Product datasheet for **GA102314**

Human IDH3G activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IDH3G CRISPRa kit - CRISPR gene activation of human isocitrate dehydrogenase (NAD(+)) 3 gamma
Symbol:	IDH3G
UniProt ID:	P51553
Synonyms:	H-IDHG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)