

Product datasheet for **GA102319**

Human IFI16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFI16 CRISPRa kit - CRISPR gene activation of human interferon gamma inducible protein 16
Symbol:	IFI16
UniProt ID:	Q16666
Synonyms:	IFNGIP1; PYHIN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)