

Product datasheet for **GA102322**

Human SP110 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SP110 CRISPRa kit - CRISPR gene activation of human SP110 nuclear body protein
Symbol:	SP110
UniProt ID:	Q9HB58
Synonyms:	IFI41; IFI75; IPR1; VODI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)