

Product datasheet for **GA102323**

Human IFIT2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFIT2 CRISPRa kit - CRISPR gene activation of human interferon induced protein with tetratricopeptide repeats 2
Symbol:	IFIT2
UniProt ID:	P09913
Synonyms:	cig42; G10P2; GARG-39; IFI-54; IFI-54K; IFI54; IFIT-2; ISG-54 K; ISG-54K; ISG54; P54
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)