

Product datasheet for **GA102325**

Human IFIT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFIT3 CRISPRa kit - CRISPR gene activation of human interferon induced protein with tetratricopeptide repeats 3
Symbol:	IFIT3
UniProt ID:	O14879
Synonyms:	CIG-49; GARG-49; IFI60; IFIT4; IRG2; ISG60; P60; RIG-G
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)