

Product datasheet for **GA102340**

Human IFNAR2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | IFNAR2 CRISPRa kit - CRISPR gene activation of human interferon alpha and beta receptor subunit 2 |
| Symbol: | IFNAR2 |
| UniProt ID: | P48551 |
| Synonyms: | IFN-alpha-REC; IFN-R; IFNABR; IFNARB; IMD45 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)